# OPERATING INSTRUCTIONS

## 08.01 SANDBOX

### Contents

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## Appendix 1: General Information

The pF-curves plotted below will be used to illustrate the soil physical characteristics that can be deduced from pF-curves. The example soil contains three different soil horizons (each of which has a known pF-Curve). These curves are referred to in Table 3.


Table 3: Determining soil characteristics from pF-Curves

Physical Characteristic	Definition and how to determine
Moisture content	Volume fraction of water filled pores at a certain matric potential. For example, at a matric potential of 1000 hPa (1 bar, pF 3.0), the A horizon has a volumetric moisture content of 20 %.
Field capacity (FC)	Moisture content at pF 2 The A horizon has a moisture content of 35 % at FC and the C horizon 24 %.
Permanent wilting point (PWP)	Moisture content at pF 4.2 The A horizon has a moisture content of 8 % at PWP and the C horizon 4 %.
Porosity	In a sandy soil, all pores are filled with water at saturation (pF0), and empty when oven-dry (pF 7). Therefore, the volume percentage between pF 0 and pF 7 is equal to the porosity in a sandy soil. In a clay soil porosity, or total pore volume, depends upon moisture content, due to swelling and shrinking. Therefore, for clay soil porosity cannot be determined from the pF -curve. The example soil is a loamy sand soil, and allows estimating porosity: At saturation, the A horizon has a volumetric moisture content of 50 %, when the soil is oven dry the moisture

	<p>content is 0 %, therefore, 50 % of the soil volume is pore space, filled with water and air, and porosity is 50 %.</p>
Volume fraction solid matter	<p>Total volume fraction minus porosity.</p> <p>Since porosity of the A horizon is 50 %, the volume fraction of pores is 0.5 and volume fraction of solid matter in the A horizon is <math>1 - 0.5 = 0.5</math></p>
Aeration status	<p>Volume of available air: porosity minus moisture content. Depending on crop type, a certain ratio between water and air supply is required for optimal crop growth.</p> <p>In the example soil, (A horizon) at a moisture potential of 1000 hPa (pF 3), moisture content is 20%, total pore space is 50 %, so the volume of available air is 30 %.</p>
Pore size distribution	<p>Shape of pF curve: Pores of similar size will be emptied at the same matric potential. The more homogenous the pore size distribution, the faster the drop in soil moisture content upon a small decrease in matric potential, and the flatter the slope of the pF curve. The steeper the slope, the more gradual the emptying of soil pores, the more heterogeneous the pore size distribution.</p> <p>In general, a heterogeneous pore size distribution is preferable for agricultural applications, since these soils have a higher water holding capacity.</p> <p>The example soil illustrates the effect of organic matter presence and biological activity in the A horizon. In the A horizon, the slope of the pF-curve is more gradual than in the C horizon, meaning that pores are emptied more gradually in the A horizon, corresponding to a heterogeneous pore size distribution. The C horizon contains a relatively large amount of pores of similar size, which are all drained around a matric potential of - 100 hPa (pF 2). A slight increase in the suction will lead to a change in moisture content of almost 10 %.</p>
Capillary conductivity	<p>The rate of capillary conductivity depends upon the amount and size of water filled pores involved in water flow. This depends upon the moisture potential of the soil.</p> <p>A decrease of the water potential (an increase in suction level) corresponds with a decrease in moisture content.</p>

	<p>Because water is forced to flow through narrow pores with a high friction, this consequently leads to a reduction in the capillary motion.</p> <p>Permeability rate depends on the distribution and amount of macro-pores.</p>
Storage capacity	<p>Storage capacity of a soil at a specific ground water level corresponds to the air volume present. Storage capability is expressed in mm water per decimetre of soil (1 mm water per 10 cm <sup>o</sup>1 volume percent).</p> <p>For the example soil, the storage capacity of the C horizon at a moisture tension of 100 hPa (pF2) is calculated as total pore space (40 %) - moisture content (25 %) = volume of air (15 %). A volumetric air content of 15 % corresponds to a storage capacity of 15 mm of water per decimetre of C horizon.</p>
Plant available soil water	<p>The amount of water between FC and PWP in volume percentage. This value should be used with caution. First, plants will start wilting with subsequent yield losses well before the permanent wilting point. Secondly, plant available soil water is replenished by capillary rise, rainfall and irrigation water.</p> <p>Eg: A fine sandy soil, rich in loam has a rooting depth of about 40 cm.</p> <ul style="list-style-type: none"> <li>• The A horizon has a depth of 20 cm.</li> <li>• The B horizon has a depth of 30 cm.</li> </ul> <p>Calculation of the amount of plant available soil water:  At field capacity, pF2.0, the A horizon will contain 35 volume % of water. At the permanent wilting point, pF4.2, the A horizon will contain 8 volume % of water.  As 1 volume % corresponds to 1 mm water per 10 cm of soil, the amount of available soil water in the A horizon is calculated as the volume % of water multiplied with the rooted depth of the soil horizon:  A horizon: 35 - 8 = 27 volume % water x 20 cm soil depth 27 x 2 dm soil depth = 54 mm  For the B horizon the calculation is similar. Notice that rooting depth is 40 cm, so roots will be present only in the upper 20 cm of the B horizon. At field capacity 27 %</p>

of water will be available, at the permanent wilting point only 6%.

B horizon:  $27 - 6 = 21$  volume % of water \* 20 cm  
rooted soil depth °  $21 * 2 = 42$  mm water

In total,  $54 + 42 = 96$  mm of water is available to plant  
growth in this particular soil.

## Appendix 2: Description of different pF-sets

To determine the soil moisture retention characteristic, the desired pF-set(s) is/are required. A balance with an accuracy of 0.01g, and a ventilated electrical drying oven (105 °C), are also necessary. Eijkelkamp supplies the following:

A universal drying-oven with temperature range 30 -220 °C, 220 V- 50 Hz.

A Sandbox for pF determination (pF0 -2.0). The Standard set for about 40 samples includes:

- Sandbox
- Containers with sand, particle size  $\pm$  73 mm, 12.5 kg each
- Filter cloth, 140 - 150 mm
- Set of 65 o-rings, diameter 49 mm: suitable for 5 cm diameter core rings

A Sand / kaolin box for pF determination (pF2.0 - 2.7). The Standard set for about 40 samples includes:

- Sand / kaolin box
- Vacuum pump and automatic suction level control system, 220 V
- Containers with synthetic sand, particle size  $\pm$  73 mm, 12.5 kg each
- Filter cloth, 140 - 150 mm
- Kaolin clay, container 2.5 kg
- Set of 65 o-rings for 5 cm diameter core rings

Pressure Membrane Apparatus (pF3.0 – 4.2). The standard set for about 15 samples includes:

- Pressure membrane extractor
- Cellophane membrane
- Soil sample retaining rings
- Filter cloth 140 - 150 mm
- Compressor 20 bar
- Air filter with support and hose

### Appendix 3: Conversion factors

100 hPa = 100 cm pressure head

= 100 cm water column

= 0.1 bar

= 0.01 Pa

= 0.01 N/m<sup>2</sup>

= 1.45 PSI

= pF (<sup>10</sup>log100)=2.0

<b>pF value</b>	<b>Matric Potential in hPa</b>	<b>Pressure in bar</b>
0	1	- 0.001
0.4	2.5	- 0.0025
1.0	10	- 0.01
1.5	31.6	- 0.0316
1.8	63.1	- 0.0631
2.0	100	- 0.1
2.3	200	- 0.2
2.7	500	- 0.5

## Appendix 4: Soil sampling

To determine the moisture retention characteristic or the pF-curve of a specific soil, *undisturbed* core samples must be collected. This is because of the major influences of both pore size distribution and soil structure on moisture retention, especially at the high matrix potentials of the operating range of suction tables.

There is no explicit prescription in literature for recommended sample sizes. Optimal sizes for core rings are determined by the size of structural elements in the soil. To obtain representative data, sample sizes should be large with respect to the size of soil aggregates, cracks, root channels or animal holes. From a practical point of view, sample diameters should not be too large as not to reduce the amount of simultaneously analysable samples, and sample height should be constrained to several centimeters; so that equilibrium conditions are reached in a reasonable period of time. According to the Dutch NEN 5787 standard, samples with a volume between 100 and 300 cm<sup>3</sup> are usually used for the suction tables, while samples with a height of more than 5 cm are discouraged, because the time needed to establish equilibrium will be long, and the accuracy of determination of pF-values near saturation will be low. In the procedures for soil analyses of the International Soil Reference and Information Centre (ISRIC), sample rings with a diameter of 5 cm and a volume of 100 cm<sup>3</sup> are recommended, while in other publications heights of 2 or 3 cm are preferred.

Eijkelpark Agrisearch Equipment recommends the use of a 100 cm<sup>3</sup> volume core ring, with an inner diameter of 50 mm (outer diameter 53 mm) and a height of 51 mm.

When pressing the core rings into the soil, care should be taken not to disturb the original setting of the soil and to completely fill the ring. Sampling conditions are best when the soil is approximately at field capacity. Ring holders may be used to facilitate insertion, especially in the subsoil. After insertion to the desired depth, the rings are carefully dug out (e.g. using the spatula provided with the Eijkelpark sample ring set), at some centimeters below the ring itself. The surplus of soil is reduced to a few millimeters, trimming it carefully with a fine iron saw, and the caps are placed on the ring for protection and to minimise evaporation losses. The remaining surplus of soil will protect the sample during transport and will be removed in the laboratory, prior to analysis. Transport the core rings in a protective case.

Since soil structure and pore size distribution have significant influence on soil water retention, several replicate samples are needed to obtain a representative pF-value. Depending on natural variability of the study area, three to six replicate samples per unit are advised.

In case the samples cannot be analysed on short notice, store the samples in a refrigerator to reduce microbial activity which might cause non-representative changes in soil structure.

- **Do not freeze the samples because soil structure will be influenced.**

## References and literature

- Klute, A. Water Retention: Laboratory Methods. IN: Methods of Soil Analysis. Part 1. Physical and Mineralogical Methods. 1986.
- Koorevaar, P., G. Menelik and C. Dirksen.  
Elements of Soil Physics. Developments in Soil Science 13/1983. Reeve, M.J. and A.D. Carter. Water Release Characteristic. IN: Soil Analysis. Physical Methods. K.A. Smith and C.E. Mullins (eds.) 1991.
- Van Reeuwijk, L.P. (ed.) Procedures for soil analysis. 1995. ISRIC Wageningen.
- Stakman, W.P., G.A. Valk and G.G. van der Harst.  
Determination of soil moisture retention curves I. 1969. ICW Wageningen. Stolte (ed.)
- Manual for soil physical measurements. Version 3. Technical Document 37. SC-DLO. 1997. Topp, G.C. and W. Zebchuk.  
The determination of soil water desorption curves for soil cores. 1979. Canadian  
Journal of Soil Science 59: 19-26.